

Комбинированные шкафы для безопасного хранения UB-T-90E PROFESSIONAL, SLIM, MODULAR

Технические характеристики

По вопросам продаж и поддержки обращайтесь:

Алматы (727)345-47-04
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Саранск (8342)22-96-24
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +7(727)345-47-04

Беларусь +(375)257-127-884

Узбекистан +998(71)205-18-59

Киргизия +996(312)96-26-47

эл.почта: auj@nt-rt.ru || сайт: <https://asecos.nt-rt.ru/>

UB-LINE | UB-T-90E PROFESSIONAL



Model UB90.060.110.2T:0014
Order No. 36521-882-38347



Model UB90.060.059.T:0014
Order No. 36511-882-38346



Model UB90.080.059.060.T:0014
Order No. 36539-882-38348

Technical data		UB90.060.110.2T:0014	UB90.060.059.T:0014	UB90.080.059.060.T:0014
External dimensions W x D x H	mm	1102 x 574 x 600	593 x 574 x 600	593 x 600 x 800
Internal dimensions W x D x H	mm	979 x 452 x 502	470 x 452 x 502	470 x 482 x 702
Weight without interior equip.	kg	180	120	143
Extraction air	DN	50	50	50
Air change 10 times	m ³ /h	2	1	2
Diff. pressure (air change 10 times EK5/AK4)	Pa	1	1	1

Main features



Incl. approved fire protection bushes



Built-in signal box for monitoring fill levels and exhaust air



Pull-out tray with positioning frame for the disposal canisters

You will find technical drawings for these cabinet models from page 228

Model	Body colour	Door colour	Order No.
UB90.060.110.2T:0014	light grey RAL 7035	stainless steel	36521-882
		light grey RAL 7035	36521-889
UB90.060.059.T:0014	light grey RAL 7035	stainless steel	36511-882
		light grey RAL 7035	36511-889
UB90.080.059.060.T:0014	light grey RAL 7035	stainless steel	36539-882
		light grey RAL 7035	36539-889

Interior equipment packages (can only be ordered in combination with the cabinet)	Material	Order No.
UB90.060.110.2T:0014		
PREMIUM		
2 x drawer	stainless steel, PE conductive	38347
1 x disposal system (3 x disposal channel)		
4 x canister 10 L		
1 x exhaust air monitoring		
BASIC		
2 x drawer	sheet steel powder-coated (electrically conductive), PE conductive	38343
1 x disposal system (3 x disposal channel)		
4 x canister 10 L		
UB90.060.059.T:0014		
PREMIUM		
1 x drawer	stainless steel, PE conductive	38346
1 x disposal system (1 x disposal channel)		
2 x canister 10 L		
1 x exhaust air monitoring		
BASIC		
1 x drawer	sheet steel powder-coated (electrically conductive), PE conductive	38341
1 x disposal system (1 x disposal channel)		
2 x canister 10 L		
UB90.080.059.060.T:0014		
PREMIUM		
1 x drawer	stainless steel, PE conductive	38348
1 x disposal system (1 x disposal channel)		
2 x canister 20 L		
1 x exhaust air monitoring		
BASIC		
1 x drawer	sheet steel powder-coated (electrically conductive), PE conductive	38345
1 x disposal system (1 x disposal channel)		
2 x canister 20 L		

Accessories for these cabinet models:

Interior equipment	Material	Order No.
UB90.060.110.2T:0014		
Castors with plinth	powder-coated sheet steel RAL 7015	31358
UB90.060.059.T:0014 / UB90.080.059.060.T:0014		
Castors with plinth	powder-coated sheet steel RAL 7015	31360
Accessories for both sides		
Fill-in set 1: Funnel small with cap (PE-HD electrical conductive) and table lead-through (GL45/GL25)	PE	37282
Fill-in set 2: ball funnel with cap (PE-HD electrical conductive) and table lead-through (GL45/GL25)		37283
Fill-in set 3: draining adapter for basin, PE-HD electrical conductive, G 1 1/2"	PE	36567

UB-LINE | UB-T-90E MODULAR



Model UB90.060.110.2T:0014
Order No. 36521-882-36825



Model UB90.060.059.T:0014
Order No. 36511-882-36824



Model UB90.080.059.060.T:0014
Order No. 36539-882-38530

Technical data

		UB90.060.059.T UB90.060.059.T:0014	UB90.060.110.2T UB90.060.110.2T:0014	UB90.080.059.060.T UB90.080.059.060.T:0014
External dimensions W x D x H	mm	593 x 574 x 600	1102 x 574 x 600	593 x 600 x 800
Internal dimensions W x D x H	mm	470 x 452 x 502	979 x 452 x 502	470 x 482 x 702
Weight without interior equip.	kg	120	180	143
Extraction air	DN	50	50	50
Air change 10 times	m ³ /h	1	2	2
Diff. pressure (air change 10 times EK5/AK4)	Pa	1	1	1

Main features



Optionally expandable in case of changing workplace conditions



Pull-out tray with positioning frame for the disposal canisters

You will find technical drawings for these cabinet models from page 228

Model	Body colour	Door colour	Order No.
Outer corpus without discharge capacity			
UB90.060.110.2T	light grey RAL 7035	light grey RAL 7035	30443-001
	laboratory white (sim. RAL 9016)	laboratory white (sim. RAL 9016)	30443-005
UB90.060.059.T	light grey RAL 7035	light grey RAL 7035	30481-001
	laboratory white (sim. RAL 9016)	laboratory white (sim. RAL 9016)	30481-005
UB90.080.059.060.T	light grey RAL 7035	light grey RAL 7035	35900-001
	laboratory white (sim. RAL 9016)	laboratory white (sim. RAL 9016)	35900-005
Outer corpus with discharge capacity			
UB90.060.110.2T:0014	light grey RAL 7035	light grey RAL 7035	36521-889
		stainless steel	36521-882
UB90.060.059.T:0014	light grey RAL 7035	light grey RAL 7035	36511-889
		stainless steel	36511-882
UB90.080.059.060.T:0014	light grey RAL 7035	light grey RAL 7035	36539-889
		stainless steel	36539-882
Interior equipment packages (can only be ordered in combination with the cabinet)		Material	Order No.
UB90.060.110.2T			
2x drawer max. interior height (V=15.5L), 2x positioning sheet		brushed stainless steel, stainless steel	36814
UB90.060.059.T			
1x drawer (disposal) left hinged (V=15.5L), 1x positioning sheet		brushed stainless steel, stainless steel	36813
UB90.080.059.060.T			
1x drawer (disposal) left hinged (V=30.0L), 1x positioning sheet		brushed stainless steel, stainless steel	37114
UB90.060.110.2T:0014			
2x drawer max. interior height (V=15.5L), 2x positioning sheet		brushed stainless steel, stainless steel	36825
UB90.060.059.T:0014			
1x drawer (disposal) left hinged (V=15.5L), 1x positioning sheet		brushed stainless steel, stainless steel	36824
UB90.080.059.060.T:0014			
1x drawer (disposal) left hinged (V=30.0L), 1x positioning sheet		brushed stainless steel, stainless steel	38530

Accessories for these cabinet models:

Interior equipment	Material	Order No.
UB90.060.110.2T(0014)		
Castors with plinth	powder-coated sheet steel RAL 7015	31358
Pipe lead-through enclosed, for customer assembly	powder-coated sheet steel red	28088
UB90.060.059.T(0014) / UB90.080.059.060.T(0014)		
Castors with plinth	powder-coated sheet steel RAL 7015	31360
Pipe lead-through enclosed, for customer assembly	powder-coated sheet steel red	28088

Please note:

The use of the pipe lead-through listed in the accessories is mandatory when installing hoses and cables of the on-site waste disposal systems. Only then is the tested fire resistance of the cabinets (type 90 according to EN 14470-1) guaranteed. The pipe lead-through is enclosed with the cabinet and can be mounted on the cabinet by the customer according to the manufacturer's specifications.

UB-LINE | UB-S-90E SLIM



Model UB90.060.110.2S
Order No. 30427-001-38354



Model UB90.060.059.S
Order No. 30472-001-38352



Model UB90.080.059.060.S
Order No. 32658-001-38357

Technical data		UB90.060.110.2S	UB90.060.059.S	UB90.080.110.060.2S	UB90.080.059.060.S
External dimensions W x D x H	mm	1102 x 574 x 600	593 x 574 x 600	1102 x 600 x 800	593 x 600 x 800
Internal dimensions W x D x H	mm	979 x 452 x 502	470 x 452 x 502	979 x 482 x 702	470 x 482 x 702
Weight without interior equip.	kg	180	120	233	143
Air change 10 times	m ³ /h	2	1	2	2
Diff. pressure (air change 10 times EK5/ AK4)	Pa	1	1	1	1
Max. load of drawers	kg	50	50	50	50

Main features



Canister incl. funnel made of high-quality, electrically conductive PE-HD

You will find technical drawings for these cabinet models from page 228

Model	Body colour	Colour	Order No.
UB90.060.110.2S	light grey RAL 7035	light grey RAL 7035	30427-001
	laboratory white (sim. RAL 9016)	laboratory white (sim. RAL 9016)	30427-005
UB90.060.059.S	light grey RAL 7035	light grey RAL 7035	30472-001
	laboratory white (sim. RAL 9016)	laboratory white (sim. RAL 9016)	30472-005
UB90.080.110.060.2S	light grey RAL 7035	light grey RAL 7035	32657-001
	laboratory white (sim. RAL 9016)	laboratory white (sim. RAL 9016)	32657-005
UB90.080.059.060.S	light grey RAL 7035	light grey RAL 7035	32658-001
	laboratory white (sim. RAL 9016)	laboratory white (sim. RAL 9016)	32658-005

Interior equipment packages (can only be ordered in combination with the cabinet)	Material	Order No.
UB90.060.110.2S		
2x drawer (V=15.5L), 2x funnel, 2x canister 10L	powder-coated sheet steel RAL 7035, PE conductive	38354
UB90.060.059.S		
1x drawer (V=15.5L), 1x funnel, 1x canister 10L	powder-coated sheet steel RAL 7035, PE conductive	38352
UB90.080.110.060.2S		
2x drawer (V=30.0L), 2x funnel, 2x canister 20L	powder-coated sheet steel RAL 7035, PE conductive	38359
UB90.080.059.060.S		
1x drawer (V=30.0L), 1x funnel, 1x canister 20L	powder-coated sheet steel RAL 7035, PE conductive	38357

Accessories for these cabinet models:

Interior equipment	Material	Order No.
UB90.060.110.2S / UB90.080.110.060.2S		
Castors with plinth	powder-coated sheet steel RAL 7015	31358
UB90.060.059.S / UB90.080.059.060.S		
Castors with plinth	powder-coated sheet steel RAL 7015	31360

UB-LINE SLIM

Disposal system for intermediate storage of waste liquids. SLIM safety storage cabinets consist of a safety storage cabinet with drawer(s) and disposal canister(s). The liquids are dispensed into the disposal canister via a screwed-on funnel.




UB-LINE SLIM

High-quality disposal canister with funnel made of electrically conductive PE-HD.



Safety storage cabinets only provide reliable protection ...


... for humankind, environment and property when they are maintained professionally! With the **asecos services** you are always on the safe side.




ARE YOU LOOKING FOR AN APPROPRIATE STORAGE SOLUTION FOR HAZARDOUS MATERIALS?

Our optimised online product configurator helps you, step by step, to find the suitable safety storage cabinet. Take the test: www.asecos-configurator.com


Top views




UB90.060.140(.UL).xxx
UB90.060.140.050(.UL).xxx




UB90.060.110(.UL).xxx
UB90.060.110.050(.UL).xxx
UB90.080.110.075.2T




UB90.060.089(.UL).xxx
UB90.060.089.050(.UL).xxx




UB90.060.059(.UL).xxx
UB90.060.059.050(.UL).xxx
UB90.060.059.T:0014



UB90.080.110.xxx



UB90.080.059.xxx
UB90.080.059.T:0014



UB90.060.110.2T:0014

По вопросам продаж и поддержки обращайтесь:

Алматы (727)345-47-04
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Саранск (8342)22-96-24
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +7(727)345-47-04

Беларусь +(375)257-127-884

Узбекистан +998(71)205-18-59

Киргизия +996(312)96-26-47

эл.почта: auj@nt-rt.ru || сайт: <https://asecos.nt-rt.ru/>